

An automated workflow for rapid alignment and identification of lipid biomarkers obtained from chip-based direct infusion nanoelectrospray tandem mass spectrometry

Jens Hoefkens (1), Tobias Kind (2), Kent Pinkerton (3), Oliver Fiehn (2)

(1) Genedata Inc, Waltham, MA; (2) UC Davis - Metabolomics, Davis, CA; (3) UC Davis Center for Health and the Environment, Davis, CA

Introduction

Polar lipids (glycerophospholipids, saccharolipids) are routinely analyzed with chromatographic and mass spectrometric techniques. Chip-based direct infusion nanoelectrospray tandem mass spectrometry can provide lipid profiles or lipid fingerprints with short infusion times (one to two minutes). Using high-resolution FT-ICR-MS or low-resolution linear ion trap mass spectra in positive and negative mode together with data dependent tandem mass spectrometry, polar lipids can be semi-quantified and identified. The major challenge is however not the analysis itself but automated data handling and data evaluation of the obtained results. We present efforts which include a reusable workflow environment for rapid alignment and identification of lipids from infusion tandem mass spectra. Results from an environmental tobacco smoke study and plasma and lung tissue are shown.


Methods

To 30 μ L plasma or 10 mg cryogenically ground lung samples 1.0 mL of chilled (isopropanol:acetonitrile:water (3:3:2; v/v/v; -20°C) solvent was added. The solution was shaken for 5 minutes in a vortexer and centrifuged 2.5 minutes at 13,000 rpm. After adding 10 μ L of 100 mM ammonium acetate to 90 μ L of the supernatant the lipid extracts were infused into a linear ion trap mass spectrometer coupled to an FT-ICR-MS using a Advion Nanomate chip-based infusion system. Mass spectra were acquired from 400-1000 m/z with positive ionization and negative ionization and 100,000 resolving power and 0V CID voltage. Data dependent MS2 scans were performed on the 40 highest peaks in the full scan spectrum and scanned in low-resolution mode with 30V CID voltage. The data alignment of the high-resolution and low-resolution mass spectra was performed using Genedata Expressionist Refiner MS. The statistical data evaluation was performed with Genedata Analyst using univariate and multivariate methods. Unique lipid species were identified by matching tandem mass spectra against a database of known fragment ions and validating results with LipidMaps tools and LipidMaps web services.

Results


Lipid extracts from heart, plasma and lungs from rats which were exposed to environmental tobacco smoke and clean air were analyzed. An alignment workflow for high-resolution and low-resolution mass spectra was built using Expressionist Refiner MS. It included spectrum grid alignment, spectrum averaging, intensity thresholding, peak detection, peak shaping and alignment export. The post-processing time for a sample setup of 70 infusion experiments took around 20 seconds. Additionally MS/MS spectra were extracted from data dependent scans. A matrix containing all samples and all aligned ions above a certain signal noise threshold was obtained as a result file. Identification of lipid species was performed by matching specific sn1 and sn2 fatty acid chain fragmentation from MS/MS data of the product ions. 50 lipid species were identified in plasma and 50 remained unknown. Different lipid patterns revealed alterations during exposure to tobacco smoke. This research project was supported by R01 ES013932 NIH/NIEHS.

Platform – nanoESI infusion with iontrap




Low-resolution and high-resolution data were obtained from a linear ion trap and a linear iontrap coupled to an Fourier transform (FT-ICR-MS) mass spectrometer. Infusion of full extracts from lung, plasma and heart samples was performed with a chip-based nanoelectrospray robot (Advion Nanomate). Infusion time per sample with data dependent MS/MS scans is two minutes to obtain about 100 mass spectra.

Mass spectral alignment with Refiner MS




Alignment of low-resolution and high-resolution mass spectral infusion data was implemented within Expressionist Refiner MS. The workflow performs peak detection, noise analysis, peak integration, isotope grouping and alignment of large datasets. The workflow runtime for 70 infusion samples is 20 seconds, allowing fast parameter optimizations. MS/MS data extraction is performed to obtain tandem mass spectra for assignment of phospholipid structures including head group and neutral loss analysis.

Alignment of infusion mass spectra




Zoom into phospholipid data region (above). Peak markers assign the correctly detected peaks. Single peak zoom before data analysis (above right) and after peak detection and alignment (lower right). Mass drifts are corrected to ensure correct assignment.

Parameter optimization and peak refinement



Depending on the mass spectral data acquisition mode, certain workflow parameters can be optimized to avoid peak detection errors or to adjust sensitivity if needed. The table shows a very high number of peaks after the first initial step. The recursive consistency filter and the signal to noise filter ensure to retain peaks which have a very low abundance but good signal-to-noise (S/N) ratio.

Data dependent MS/MS for class analysis



The analysis of the different phospholipid classes including Glycerophosphocholines (GPCho), Sphingomyelins (SM), Glycerophosphoethanolamines (GPEtn), Glycerophosphoinositols (GPIins), Glycerophosphoserines (GPSer) and others can be performed by subtraction of specific neutral losses from the precursor ion in positive and negative mode. The precursor mass is then used to compute or lookup only species from a specific lipid class. In case of overlaps or unclear spectra multiple species can occur.

sn1 and sn2 fatty acid analysis from MS/MS

